

RPI Part #	Industry Part #	Description	Section 8 Page #
<b>LAMPS &amp; BULBS</b>			
AMB028	NE51	LIGHT BULB	201
AIB041	GE345	BULB, NEON	201
LMP001	GE252	LAMP	200
LMP002	GE778	LAMP	200
LMP003	GE787	LAMP	200
LMP004	GE1460	LAMP	200
LMP005	15S11/102	LAMP	200
LMP006	FCR7023	LAMP	200
LMP007	ESA/FHD	LAMP	200
LMP008	ESB7388	LAMP	200
LMP009	GE1493	LAMP	200
LMP010	15S11/CL13	LAMP	200
LMP011	BRL7027	LAMP	200
LMP012	13347W	LAMP	200
LMP013	GE1634	LAMP	200
LMP014	EXN	LAMP	200
LMP015	25T10	LAMP	200
LMP101	GE1004	SPECTROPHOTOMETER LAMP	200
LMP102	D2A	SPECTROPHOTOMETER LAMP	200
LMP103	GE82	SPECTROPHOTOMETER LAMP	201
LMP104	BWD	SPECTROPHOTOMETER LAMP	201
LMP600	FCS	LAMP	201
LMP601	H-3	LAMP	201
LMP602	125T10/P	LAMP	201
LMP603	EJL	LAMP	201
LMP604	ENZ	LAMP	201
LMP605	ETJ	LAMP	201
LMP606	ELZ	LAMP	201
LMP607	FDS/DZE	LAMP	201
LMP608	75T3Q/CL	LAMP	201
PCL606	150T4Q/CL	LAMP	201
PCT685	N/A	SOCKET RETRACTOR TOOL	201

● NEW PART

Call Us Toll-Free Monday–Friday, 8:00 AM–4:30 PM (Pacific Time), (800) 221-9723

SECTION

8

PAGE

100

UPDATED

6/98

Lamps & Bulbs

NOTES:

Lined area for handwritten notes.

**RPI Part #LMP001**  
**Industry Part #GE252**

- LAMP**
- 2.5V, .35A
  - Base: Single contact, midget flanged




**RPI Part #LMP007**  
**Industry Part #ESA/FHD**

- LAMP**
- 6V, 10W
  - Base: 2-pin




**RPI Part #LMP013**  
**Industry Part #GE1634**

- LAMP**
- 20V, 1.0A
  - Base: Double contact, prefocus
  - Halogen



**RPI Part #LMP002**  
**Industry Part #GE778**

- LAMP**
- 6.0V, 20W
  - Base: 2-pin
  - Halogen




**RPI Part #LMP008**  
**Industry Part #ESB 7388**

- LAMP**
- 6V, 20W
  - Base: 2-pin



**RPI Part #LMP014**  
**Industry Part #EXN**

- LAMP**
- 12V, 50W
  - Base: 2-pin
  - Incandescent




**RPI Part #LMP003**  
**Industry Part #GE787**

- LAMP**
- 6.0V, 10W
  - Base: 2-pin
  - Halogen




**RPI Part #LMP009**  
**Industry Part #GE1493**

- LAMP**
- 6.5V, 2.75A
  - Base: Double contact, bayonet
  - Halogen




**New RPI Part #LMP015**  
**Old RPI Part #HAL007**  
**Industry Part #25T10**

- LAMP**
- 120V, 25W
  - Base: Medium screw, prefocus



**RPI Part #LMP004**  
**Industry Part #GE1460**

- LAMP**
- 6.5V, 2.75A
  - Base: Double contact, prefocus
  - Halogen




**RPI Part #LMP010**  
**Industry Part #15S11/CL13**

- LAMP**
- 115V, 15W
  - Base: Candelabra




**RPI Part #LMP101**  
**Industry Part #GE1004**

- SPECTROPHOTOMETER LAMP**
- 12.8V, .94A
  - Base: Double contact, bayonet with mounting bracket
  - Halogen
  - 18 GA black and red wire with connectors



**RPI Part #LMP005**  
**Industry Part #15S11/102**

- LAMP**
- 120V, 15W
  - Base: Medium screw, prefocus




**RPI Part #LMP011**  
**Industry Part #BRL 7027**

- LAMP**
- 12.V, 50W
  - Base: Double contact, miniature



**RPI Part #LMP102**  
**Industry Part #D2A**

- SPECTROPHOTOMETER LAMP**
- 10.8V, 30W
  - Base: Mini, 2-pin



**RPI Part #LMP006**  
**Industry Part #FCR 7023**

- LAMP**
- 12V, 100W
  - Base: 2-pin
  - Halogen



**RPI Part #LMP012**  
**Industry Part #13347W**

- LAMP**
- 6V, 15W
  - Base: Double contact, bayonet candelabra
  - Halogen



Lamps & Bulbs

RPI Part #LMP103  
Industry Part #GE82

**SPECTROPHOTOMETER LAMP**  
• 6.5V, 1.02A  
• Base: Double contact, bayonet  
• Halogen




RPI Part #LMP104  
Industry Part #BWD

**SPECTROPHOTOMETER LAMP**  
• 8.5V, 4.2A  
• Base: Double contact, bayonet with mounting bracket  
• 18 GA black and red wire 8" with connectors




RPI Part #LMP600  
Industry Part #FCS

**LAMP**  
• 24V, 150W  
• Base: 2-pin




RPI Part #LMP601  
Industry Part #H-3

**LAMP**  
• 12V, 55W  
• Base: Wire  
• Halogen




RPI Part #LMP602  
Industry Part #125T10P

**LAMP**  
• 120V, 125W  
• Base: Medium, prefocus




RPI Part #LMP603  
Industry Part #EJL

**LAMP**  
• 24V, 200W  
• Base: 2-pin




RPI Part #LMP604  
Industry Part #ENZ

**LAMP**  
• 30V, 50W  
• Base: 2-pin




RPI Part #LMP605  
Industry Part #ETJ

**LAMP**  
• 120V, 250W  
• Base: 2-pin, special




RPI Part #LMP606  
Industry Part #ELZ

**LAMP**  
• 21V, 150W  
• Base: 2-pin




RPI Part #LMP607  
Industry Part #FDS/DZE

**LAMP**  
• 24V, 150W  
• Base: 2-pin




RPI Part #LMP608  
Industry Part #75T3Q/CL

**LAMP**  
• 28V, 75W  
• Base: Double contact, bayonet  
• Halogen



RPI Part #AMB028  
Industry Part #NE51

**LIGHT BULB**  
• Standard NE51 bulb




RPI Part #AIB041  
Industry Part #GE345

**BULB, NEON**  
• T3, 1/4 E10 base  
• 2 watts, 24/30 VDC




RPI Part #PCL606  
Industry Part #150T4Q/CL

**LAMP**  
• 24V, 150W  
• Base: Double ended  
• Halogen



RPI Part #PCT685  
**SOCKET RETRACTOR TOOL**

• For installation of lamp (RPI Part #PCL606)



Use the handy fax order form located at the front of the catalog.