


Models: DV5, DV7, DV10, DV12, SC3, SC5, SC7, SC10, & SC12

New PM Kits to fit

**A-dec®/Dean®
Dental Compressors
and Dental Vacuum Units!**

PM Kits to fit A-dec/Dean Dental Compressors




RPI Part #CMK293
OEM Part #A0050

COMPRESSOR PM KIT

- For Oil-less dual, triple and quad headed compressors
- All parts also sold separately
- Includes all parts as shown plus a PM Sticker and PM checklist
- See *the RPI Advantage* below

Fits: Compressor Body
Models: SC3, SC5 & SC7



RPI Part #CMK294
OEM Part #A0150

COMPRESSOR PM KIT

- For Oil-less 5 and 6 headed compressors
- All parts also sold separately
- Includes all parts as shown plus a PM Sticker and PM checklist
- See *the RPI Advantage* below

Fits: Compressor Body
Models: SC10 & SC12

PM Kit to fit A-dec/Dean Dental Vacuum Units




RPI Part #VPK183
OEM Part #A0049

VACUUM PM KIT

- **Service Tip:** Air filter life depends on vacuum usage, nature of environment, and amount of aspirated material. Inspect air filter every 1,000 hours and replace every 2,000 hours or when visibly dirty.
- All parts also sold separately
- Includes all parts as shown plus a PM Sticker and PM checklist
- See *the RPI Advantage* below

Fits: Inside Separator Tank and Drain Assembly
Models: DV5, DV7, DV10 & DV12




RPI Filters and PM Kits include a PM Sticker that can be attached to the compressor or vacuum unit for recording the date that the filters were replaced. It's the RPI Advantage!

rpi replacement parts industries, inc.
CALL (800) 221-9723 • FAX (818) 882-7028 • EMAIL order@riparts.com • WEBSITE www.riparts.com

The above parts are manufactured by Replacement Parts Industries, Inc. to fit the equipment noted. All product names used in this document are trademarks or registered trademarks of their respective holders. 2010AD Show 02/21

Air Filters to fit A-dec/Dean Dental Compressors!



RPI Part #CMF295


OEM Part # (OEM Part # Not Available)

AIR FILTER

- Size: 9.87" lg. x 3.25" wide
- Material: 60 PPI Quadrafoam, color black, aluminum frame and stainless steel mesh
- Included in: Compressor PM Kit (RPI Part #CMK293) and Compressor PM Kit (RPI Part #CMK294)
- See *the RPI Advantage* on reverse side

Fits: Compressor Body

Models: SC3, SC5, SC7, SC10 & SC12



RPI Part #CMF296

OEM Part # (OEM Part # Not Available)

AIR FILTER

- Size: 18.44" lg. x 5.44" wide
- Material: 25 PPI Quadrafoam, color green, aluminum frame and stainless steel mesh
- Included in: Compressor PM Kit (RPI Part #CMK293) and Compressor PM Kit (RPI Part #CMK294)
- See *the RPI Advantage* on reverse side

Fits: Compressor Body

Models: SC3, SC5, SC7, SC10 & SC12


New Parts to fit A-dec/Dean Dental Compressors & Vacuums!

ADDRESS SERVICES REQUESTED

PERMIT NO. 2399
VAN NUYS, CA
U.S. POSTAGE PAID
PRESORTED STANDARD

rpi!
replacement parts industries, inc.
625 Cochran Street • Simi Valley, California • 93065

Parts to fit A-dec/Dean Dental Vacuum Units!



RPI Part #VPE184


OEM Part # (OEM Part # Not Available)

FILTER ELEMENT

- 5 micron
- Size 7.88" OD x 4.81" lg.
- Included in: Vacuum PM Kit (RPI Part #VPK183)
- See *the RPI Advantage* on reverse side

Fits: Inside Separator Tank

Models: DV5, DV7, DV10 & DV12



RPI Part #VPV120

OEM Part # (OEM Part # Not Available)

CHECK VALVE (1-1/2")

- Ports: 1-1/2" FPT
- Material: Clear PVC
- Size: 3.38" dia. x 5.75" lg.
- Included in: Vacuum PM Kit (RPI Part #VPK183)

Fits: Inside Separator Tank and Drain Assembly

Models: DV5, DV7, DV10 & DV12