



DentalEZ[®] Osprey[™] Dental Compressors!

Models: Osprey 13, Osprey 22, Osprey 23, Osprey 24, Osprey 25 & Osprey 28

COALESCING & PARTICULATE FILTER ASSEMBLIES




RPI Part #CMA229
OEM Part #004960

COALESCING FILTER ASSEMBLY

- Includes all parts as shown
- Also available: Mounting Bracket (RPI Part #CMB262)

Fits: Dryer Assembly
Models: Osprey 13, Osprey 22, Osprey 23, Osprey 24, Osprey 25 & Osprey 28




RPI Part #CMA230
OEM Part #004961

PARTICULATE FILTER ASSEMBLY

- Includes all parts as shown
- Also available: Mounting Bracket (RPI Part #CMB262)

Fits: Dryer Assembly
Models: Osprey 13, Osprey 22, Osprey 23, Osprey 24, Osprey 25 & Osprey 28



The RPI Advantage - To help make your job easier when replacing the Coalescing and Particulate Filter Assemblies, RPI is proud to sell the individual components separately, as well as in Kits. It's the RPI Advantage!




RPI Part #CME227
OEM Part #005150

COALESCING FILTER ELEMENT

- 0.30 Micron filter material
- Includes all parts as shown
- It is recommended that this filter be replaced annually.

Fits: Coalescing Filter Assembly
Models: Osprey 13, Osprey 22, Osprey 23, Osprey 24, Osprey 25 & Osprey 28




RPI Part #CME228
OEM Part #005151

PARTICULATE FILTER ELEMENT

- 0.01 Micron filter material
- Includes all parts as shown
- It is recommended that this filter be replaced annually.

Fits: Particulate Filter Assembly
Models: Osprey 13, Osprey 22, Osprey 23, Osprey 24, Osprey 25 & Osprey 28




RPI Part #CMB232
OEM Part # (No OEM Part # Available)

BOWL (AUTO DRAIN)

- Includes: Float assembly, auto drain, and (1 pc) Bowl Gasket (RPI Part #CMG254)
- See *The RPI Advantage* above

Fits: Coalescing Filter Assembly
Models: Osprey 13, Osprey 22, Osprey 23, Osprey 24, Osprey 25 & Osprey 28




RPI Part #CMB233
OEM Part # (No OEM Part # Available)

BOWL (MANUAL DRAIN)

- Includes: Float assembly, manual drain, and (1 pc) Bowl Gasket (RPI Part #CMG254)
- See *The RPI Advantage* above

Fits: Particulate Filter Assembly
Models: Osprey 13, Osprey 22, Osprey 23, Osprey 24, Osprey 25 & Osprey 28




RPI Part #CMI231
OEM Part # (No OEM Part # Available)

DIFFERENTIAL PRESSURE INDICATOR

- Includes: Gasket and mounting screws
- See *The RPI Advantage* above

Fits: Coalescing Filter Assembly and Particulate Filter Assembly
Models: Osprey 13, Osprey 22, Osprey 23, Osprey 24, Osprey 25 & Osprey 28



RPI Part #CMG254
OEM Part # (No OEM Part # Available)

BOWL GASKET

- 3 per package
- Size: 2.200" ID x 2.400" OD x .100" SQ. C/S
- See *The RPI Advantage* above

Fits: Dryer Assembly, Coalescing Filter Assembly and Particulate Filter Assembly
Models: Osprey 13, Osprey 22, Osprey 23, Osprey 24, Osprey 25 & Osprey 28



CMG254
CROSS SECTION
SHOWN ACTUAL SIZE




RPI Part #CMB262
OEM Part # (No OEM Part # Available)

MOUNTING BRACKET

- Includes: Mounting screws (M5 x 10)
- Use this bracket to externally mount the Coalescing Filter Assembly and Particulate Filter Assembly to the wall. It's the RPI Advantage!

Models: Osprey 13, Osprey 22, Osprey 23, Osprey 24, Osprey 25 & Osprey 28



replacement parts industries, inc.

CALL (800) 221-9723 • FAX (818) 882-7028 • EMAIL order@ripiparts.com • WEBSITE www.ripiparts.com

The above parts are manufactured by Replacement Parts Industries, Inc. to fit the equipment noted. All product names used in this document are trademarks or registered trademarks of their respective holders. Show 08/18



COMPRESSOR PM KITS

It is recommended that these PM Kits be replaced annually.


RPI Part #CMK223
OEM Part #005147SP



COMPRESSOR PM KIT

- For 2 air-intake filter compressors
- Includes all parts as shown

Models: Osprey 13, Osprey 22, & Osprey 23




RPI Part #CMK224
OEM Part #005148SP



COMPRESSOR PM KIT

- For 4 air-intake filter compressors
- Includes all parts as shown

Model: Osprey 25




RPI Part #CMK225
OEM Part #005149SP



COMPRESSOR PM KIT

- For 6 air-intake filter compressors
- Includes all parts as shown

Model: Osprey 28




RPI Part #CMK226
OEM Part #005350SP



COMPRESSOR PM KIT

- For 3 air-intake filter compressors
- Includes all parts as shown

Model: Osprey 24



DRYERS

It is recommended that these Dryers be replaced every 5 years.

RPI Part #CMD234
OEM Part #004018SP




DRYER

- 2 per package
- Size: 12.34" lg. x 2.00" OD
- Includes all parts as shown
- **Notes:** There are two dryers per assembly.

Fits: Dryer Assembly

Models: Osprey 13, Osprey 22 & Osprey 23



RPI Part #CMD235
OEM Part #004020SP




DRYER

- 2 per package
- Size: 14.25" lg. x 2.52" OD
- Includes all parts as shown
- **Notes:** There are two dryers per assembly.

Fits: Dryer Assembly

Models: Osprey 24 & Osprey 25



RPI Part #CMD236
OEM Part #005413SP




DRYER

- 2 per package
- Size: 15.50" lg. x 3.00" OD
- Includes all parts as shown
- **Notes:** There are two dryers per assembly.

Fits: Dryer Assembly

Model: Osprey 28



RPI Part #RX0100
OEM Part #
(No OEM Part # Available)




O-RING

- 6 per package
- Size: .549" ID x .755" OD x .103" C/S

Fits: Dryer Assembly

Models: Osprey 13, Osprey 22, Osprey 23, Osprey 24, Osprey 25 & Osprey 28



INLET & OUTLET SHUTTLE DISKS

RPI Part #CMD239

OEM Part #


(No OEM Part # Available)

OUTLET SHUTTLE DISK

- Material: Teflon®
- Size: .57" OD dia x .25" thk

Fits: Dryer Assembly

Models: Osprey 13, Osprey 22 & Osprey 23



RPI Part #CMD240

OEM Part #


(No OEM Part # Available)

OUTLET SHUTTLE DISK

- Material: Teflon®
- Size: .57" OD x .25" thk

Fits: Dryer Assembly

Models: Osprey 24 & Osprey 25



RPI Part #CMD241

OEM Part #


(No OEM Part # Available)

OUTLET SHUTTLE DISK

- Material: Teflon®
- Size: .57" OD x .25" thk

Fits: Dryer Assembly

Models: Osprey 28



RPI Part #CMD242

OEM Part #


(No OEM Part # Available)

INLET SHUTTLE DISK

- Material: Teflon®
- Size: .57" OD x .25" thk


Fits: Dryer Assembly

Models: Osprey 13, Osprey 22, Osprey 23, Osprey 24, Osprey 25 & Osprey 28



SOLENOID VALVES & REPAIR KITS

RPI Part #CMV255
OEM Part #004980



SOLENOID VALVE


- Ports: 1/8" FPT; 2-Way; N.O.
- 24VDC; 4 1/2 W; (128 Ohms)
- Also available: Purge Muffler (RPI Part #CMM261) and Wire Connector with Gasket (RPI Part #TUC084)

Fits: Dryer Assembly

Models: Osprey 13, Osprey 22, Osprey 23, Osprey 24, Osprey 25 & Osprey 28

PURGE MUFFLERS

RPI Part #CMM261
OEM Part #004856



PURGE MUFFLER

- Port Size: 1/8" NPT
- Size: .51" OD x 1.50" lg.
- Material: Sintered copper alloy

Fits: Dryer Assembly


Models: Osprey 13, Osprey 22, Osprey 23, Osprey 24, Osprey 25 & Osprey 28

SOLENOID VALVES & REPAIR KITS

It's the **rpi** Advantage!

When servicing the solenoid valve, you now have a choice to replace the entire valve or just repair it with RPI's Solenoid Valve Repair Kit. It's the RPI Advantage!

RPI Part #CMV257
OEM Part #8801-024




SOLENOID VALVE

- Ports: 1/4" FPT; 2-Way; N.C.
- 24VAC 50/60 Hz., 8W (10 Ohms)
- **Note:** There are 2 solenoid valves per unit
- Also available: Solenoid Valve Repair Kit (RPI Part #CMK258), Purge Muffler (RPI Part #CMM237) and Wire Connector with Gasket (RPI Part #TUC084)

Fits: Dryer Assembly

Models: Osprey 13, Osprey 22, Osprey 23, Osprey 24 & Osprey 25

RPI Part #CMM237
OEM Part #004088



PURGE MUFFLER

- Port Size: 1/4" MPT
- Size: 13/16" OD x 2-7/16" lg
- Material: Polyvinyl foam

Fits: Dryer Assembly

Models: Osprey 13, Osprey 22, Osprey 23, Osprey 24 & Osprey 25


RPI Part #CMK258
OEM Part #
(No OEM Part # Available)

SOLENOID VALVE REPAIR KIT


- Buna-N Seal
- Includes all parts as shown
- **Note:** There are 2 solenoid valves per unit
- Also available: Solenoid Valve (RPI Part #CMV257)

Fits: Dryer Assembly

Models: Osprey 13, Osprey 22, Osprey 23, Osprey 24 & Osprey 25



RPI Part #CMV259
OEM Part #8801-025




SOLENOID VALVE

- Ports: 3/8" FPT; 2-Way; N.C.
- 24VDC 50/60 Hz., 11.5W (50 Ohms)
- **Note:** There are 2 solenoid valves per unit
- Also available: Solenoid Valve Repair Kit (RPI Part #CMK260), Purge Muffler (RPI Part #CMM238) and Wire Connector with Gasket (RPI Part #TUC084)

Fits: Dryer Assembly

Models: Osprey 28

RPI Part #CMM238
OEM Part #005082



PURGE MUFFLER

- Port Size: 3/8" NPT
- Size: 1" OD x 2 5/8" lg
- Material: Polyvinyl foam

Fits: Dryer Assembly

Models: Osprey 28


RPI Part #CMK260
OEM Part #
(No OEM Part # Available)

SOLENOID VALVE REPAIR KIT


- Includes all parts as shown
- **Note:** There are 2 solenoid valves per unit
- Also available: Solenoid Valve (RPI Part #CMV259)

Fits: Dryer Assembly

Models: Osprey 28



RPI Part #TUC084
OEM Part #
(No OEM Part # Available)




WIRE CONNECTOR WITH GASKET

- Includes all parts as shown

Fits: Dryer Solenoid Valves

Models: Osprey 13, Osprey 22, Osprey 23, Osprey 24, Osprey 25 & Osprey 28

RPI Part #CMV217
OEM Part #87045 & 004761SP



CHECK VALVE


- Ports: 3/8" FPT Inlet x 3/8" MPT Outlet
- Size: 7/8" Hex x 1.94" lg.

Fits: Manifold Assembly

Models: Air Techniques: AirStar10MD, AirStar12MD, AirStar21MD, AirStar22MD, AirStar30MD, AirStar40MD, AirStar50MD, AirStar70MD & AirStar100MD;

Dental-EZ: Osprey 13, Osprey 22, Osprey 23, Osprey 24, Osprey 25 & Osprey 28

RPI Part #RXF024
OEM Part #003863SP



ELBOW (3/8" PUSH-IN x 3/8" MPT)

- 2 Per package
- Material: Brass swivel fitting with plastic chuck

Fits: Dryer Assembly

Models:

Air Techniques: AS70;
ODental-EZ: Osprey 13, Osprey 22, Osprey 23, Osprey 24, Osprey 25 & Osprey 28

New Line of Parts to fit

DentalEZ[®] Osprey[™] Dental Compressors

- Bowl (Auto Drain)
- Bowl (Manual Drain)
- Check Valve
- Coalescing Filter Assembly
- Coalescing Filter Element
- Compressor PM Kits
- Differential Pressure Indicator
- Dryers
- Inlet Shuttle Disk
- Outlet Shuttle Disks
- Particulate Filter Assembly
- Particulate Filter Element
- Purge Mufflers
- Solenoid Valves
- Solenoid Valve Repair Kits
- and more !

Fits Osprey 13, Osprey 22, Osprey 23, Osprey 24, Osprey 25 & Osprey 28

New Parts to fit DentalEZ[®] Dental Compressors!

ADDRESS SERVICES REQUESTED

PERMIT NO. 2399
VAN NUYS, CA
U.S. POSTAGE PAID
PRESORTED STANDARD


replacement parts industries, inc.
625 Cochran Street • Simi Valley, California • 93065



New PM KITS to fit DentalEZ[®] Osprey[™] Dental Compressors!

Is it that time again? When was the last time you replaced the filter elements in your DentalEZ Osprey dental compressor?

Keep your equipment running in peak performance with RPI's new PM Kits. Each Kit includes the Coalescing Filter Element, Particulate Filter Element and Air Filter Elements that should be replaced on an annual basis for your model - in one convenient kit. Order one, then you're done!



Coalescing
Filter Element
(RPI Part #CME227)

Particulate
Filter Element
(RPI Part #CME228)

Air
FilterElement
(RPI Part #CME067)

Fits Osprey 13, Osprey 22, Osprey 23, Osprey 24, Osprey 25 & Osprey 28