
Dental Equipment



replacement parts industries, inc.

CALL (800) 221-9723 • FAX (818) 882-7028 • EMAIL order@rpiparts.com • WEBSITE www.rpiparts.com

The above parts are manufactured by Replacement Parts Industries, Inc. to fit the equipment noted. All product names used in this document are trademarks or registered trademarks of their respective holders. Show 06/15 Updated 07/15

SECTION

4

PAGE

701

UPDATED

06/15

ENTAL VACUUM UN

Air Techniques, Apollo/Midmark®, Dental EZ®/Custom Air/RamVac, Matrx/Midmark®, Midmark® & Tech West



RPI Part #GOG046

rpi


OEM Part # See Chart on Back Page

- VACUUM GAUGE
- Size: 1-1/2" dia.
- Dual Scale: 0 to 30 In. Hg (0 to 1 Bar)
- Material: Stainless steel case and bezel
- Port: 1/8" MPT center back mount
- Fits: Manifold Assembly



RPI Part #VPG044 OEM Part # See Chart on Back Page

- VACUUM GAUGE (LIQUID FILLED)
- Size: 1-1/2" dia.
- Dual Scale: 0 to 30 In. Hg (0 to 1 Bar)
- Material: Stainless steel case and bezel
- Port: 1/8" MPT center back mount
- Fits: Manifold Assembly



RPI Part #VPG045 OEM Part # See Chart on Back Page

VACUUM GAUGE

• Size: 2.0" dia.

• Size: 2.0" dia.


• Dual Scale: 0 to 30 In. Hg (0 to 1 Bar)

• Ports: 1/4" MPT bottom mount

Fits: Manifold Assembly

Material: Stainless steel case and bezel

- Dual Scale: 0 to 30 In. Hg (0 to 1 Bar)
- Material: Stainless steel case and bezel
- Port: 1/8" MPT center back mount
- Includes all parts as shown
- Fits: Manifold Assembly



RPI Part #VPG046 OEM Part # See Chart on Back Page


VACUUM GAUGE (LIQUID FILLED)

- Size: 2.0" dia.
- Dual Scale: 0 to 30 In. Hg (0 to 1 Bar)
- Material: Stainless steel case and bezel
- Port: 1/8" MPT center back mount
- Includes all parts as shown
- Fits: Manifold Assembly



- Size: 2.0" dia.
- Dual Scale: 0 to 30 In. Hg (0 to 1 Bar)
- Material: Stainless steel case and bezel
- Port: 1/4" MPT bottom mount
- Fits: Manifold Assembly






P

Air Techniques, Apollo/Midmark®, Dental EZ®/Custom Air/RamVac, Matrx/Midmark®, Midmark® & Tech West 0 **RPI Part #VPV058 RPI Part #VPV057 OEM Part #** See Chart on Back Page **OEM Part #** See Chart on Back Page CHECK VALVE (1/4" IN-LINE) RPI Part #VPG039 Cracking Pressure: 1.0 PSI CHECK VALVE (3/4") **OEM Part #** See Chart on Back Page • Ports: 1/4" FPT Inlet x 1/4" MPT Outlet • Ports: 3/4" Slip x 3/4" Slip Schedule 40 • Material: White PVC plastic • Material: Brass body, Buna-N (Nitrile) seal MANIFOLD GASKET Fits: Manifold Assembly Fits: Manifold Assembly Fits: Manifold Assembly **RPI Part #VPV051 OEM Part # RPI Part #VPG040** See Chart on Back Page **OEM Part #** See Chart on Back Page CHECK VALVE (3/4") MANIFOLD GASKET Ports: 3/4" FPT Fits: Manifold Assembly Material: Bronze body **RPI Part #VPV055** The RPI Check Valve is a rebuildable valve **OEM Part #** See Chart on Back Page with a Buna-N (Nitrile) seal - It's the RPI CHECK VALVE (3/8") Advantage! Ports: 3/8" FPT Also available: Check Valve Repair Kit Material: Bronze body 0 (RPI Part #VPK052) • The RPI Check Valve is a rebuildable valve Fits: Manifold Assembly with a Buna-N (Nitrile) seal - It's the RPI 0 Advantage! Also available: Check Valve Repair Kit **RPI Part #VPG041** (RPI Part #VPK050) **OEM Part #** See Chart on Back Page Fits: Manifold Assembly MANIFOLD GASKET **RPI Part #VPK052** Fits: Manifold Assembly **OEM Part #** See Chart on Back Page **CHECK VALVE REPAIR KIT (3/4")** • Material: Buna-N (Nitrile) seat Fits: Only Fits RPI's Check Valve (3/4") (RPI Part #VPV051) **RPI Part #VPG042 OEM Part #** See Chart on Back Page **RPI Part #VPV049** MANIFOLD GASKET **OEM Part #** See Chart on Back Page Fits: Manifold Assembly **RPI Part #VPV053** CHECK VALVE (1/2") **OEM Part #** • Ports: 1/2" FPT See Chart on Back Page Material: Bronze body **CHECK VALVE (1")** • The RPI Check Valve is a rebuildable valve Ports: 1" FPT with a Buna-N (Nitrile) seal - It's the RPI Material: Bronze body Advantage! • The RPI Check Valve is a rebuildable valve Also available: Check Valve Repair Kit with a Buna-N (Nitrile) seal - It's the RPI (RPI Part #VPK050) **RPI Part #VPG043** Advantage! Fits: Manifold Assembly **OEM Part #** See Chart on Back Page • Also available: Check Valve Repair Kit (RPI Part #VPK054) MANIFOLD GASKET Fits: Manifold Assembly Fits: Manifold Assembly **RPI Part #VPK050 RPI Part #VPK054 OEM Part #** See Chart on Back Page **OEM Part #** See Chart on Back Page CHECK VALVE REPAIR KIT (1/2" & 3/8") **RPI Part #VPG093** • Material: Buna-N (Nitrile) seal **CHECK VALVE REPAIR KIT (1") OEM Part #** See Chart on Back Page Fits: Only Fits Check Valve (1/2") (RPI • Material: Buna-N (Nitrile) seat Part #VPV049) and Check Valve (3/8") Fits: Only Fits RPI's Check Valve (1") MANIFOLD GASKET (RPI Part #VPV055) (RPI Part #VPV053) Fits: Manifold Assembly






- Includes: Sink Adapter (RPI Part #VPA084) and Cuspidor Drain Valve (RPI Part #VPV083)
- Also available: Diaphragm (RPI Part #VPD086) and Bowl Kit (Long) (RPI Part #VPK088)

Fits: Standard or Operatory Sinks

Note: Adapts all OEM sinks to a central vacuum line when there is not a normal gravity drain available.



Material: Silicone rubber

Valve (RPI Part #VPV083)

• Size: 2-9/16" dia. without mounting holes

Also available: Sink Drain Valve Assembly

(RPI Part #VPK082) and Cuspidor Drain

Fits: Sink and Cuspidor Drain Valve Assemblies (Vada Trap) Serial Number 0001 and

above (Including All Vada Traps Sold By RPI)



replacement parts industries, inc. 20338 Corisco Street • Chatsworth, California • 91311

ADDRESS SERVICES REQUESTED

PRESORTED STANDARD U.S. POSTAGE PAID VAN NUYS, CA PERMIT NO. 2399

List of New RPI Parts to fit Dental Vacuum Units See our website www.rpiparts.com for specific models these parts fit!

RPI		Air	Apollo/	Dental <i>EZ</i> ®/Custom				Vada
Part #	Description	Techniques	Midmark®	Air/RamVac	Midmark®	Midmark®	Tech West	Traps
CMF013	Filter Element		20207400 & 202047500					
G0G046	Vacuum Gauge	55616			62947800	62947800		
RPF936	Female Connector			No OEM Part # Available				
	(1/4" Poly x 1/4" FPT) - 2/pkg							
RP0934	0-ring - 3/pkg	No OEM Part # Available						
VPA084	Sink Adaptor	55000	D0050705	04545040	00007001	00007001	VEA 40	No OEM Part # Available
VPC069 VPC070	Solids Collector (3/4") Solids Collector (1")	55862	PCC50705 PCV50706	64545040 64501940	62937301	62937301 77005147	VFA-40	
VPC070 VPC071	Solids Collector (1-1/4")	55612	PCV50708	04301940	77005155	77005155		
VPD085	Diaphragm	55012	F 6V 307 00		11003133	77003133		No OEM Part # Available
VPD086	Diaphragm							No OEM Part # Available
VPE064	Filter Element and Screen	55713						
VPE065	Filter Element				62935604	62935604		
VPE067	Filter Element and Screen			64568145				
VPE068	Filter Element and Screen						No OEM Part # Available	
VPF066	Filter Assembly			64568159				
VPG039	Manifold Gasket	55617						
VPG040	Manifold Gasket				61203200	61203200		
VPG041	Manifold Gasket		MMS71167	04500107				
VPG042	Manifold Gasket			64568137			MG-100	
VPG043 VPG044	Manifold Gasket Vacuum Gauge (Liquid Filled)	55616			62947800	62947800	IVIG-100	
VPG044 VPG045	Vacuum Gauge (Liquid Filled)	55680			62947800	62947600		
VPG045 VPG046	Vacuum Gauge (Liquid Filled)	55680						
VPG047	Vacuum Gauge	33000	014-0550-00	64568188				
VPG048	Vacuum Gauge (Liquid Filled)		014-0550-00	64568188			VPG-100	
VPG080	Bowl Gasket - 3/pkg	55625		No OEM Part # Available	77007018	77007018	VFG-100	
VPG081	Bowl Gasket - 3/pkg	No OEM Part # Available	MRP70936	No OEM Part # Available	77005158	77005158		
VPG091	Bowl Gasket							No OEM Part # Available
VPG093	Manifold Gasket		MMS71165					
VPG094	Vacuum Gauge			64568134				
VPG095	Vacuum Gauge (Liquid Filled)			64568134				
VPK050	Check Valve Repair Kit (1/2" & 3/8")		No OEM Part # Available	No OEM Part # Available				
VPK052	Check Valve Repair Kit (3/4")	No OEM Part # Available	No OEM Part # Available	No OEM Part # Available	No OEM Part # Available	No OEM Part # Available		
VPK054	Check Valve Repair Kit (1")		No OEM Part # Available	No OEM Part # Available				
VPK062	Vacuum Relief Valve Repair Kit Solids Collector	55000	AV/AC0001	No OEM Part # Available	77005107	77005167	VFB-100	
VPK072	Replacement Kit	55880	AVA60001	No OEM Part # Available	77005167			
VPK073	Solids Collector Replacement Kit (Case of 12)	No OEM Part # Available	No OEM Part # Available	No OEM Part # Available	77005177	77005177	No OEM Part # Available	
VPK074	Solids Collector Replacement Kit	55094	AVA600012	64568122	77005180	77005180		
VPK075	Solids Collector Replacement Kit (Case of 8)	No OEM Part # Available	No OEM Part # Available	No OEM Part # Available	77005181 & 77005178	77005181 & 77005178		
VPK082	Sink Drain Valve Assembly				77003170	77003170		45R062
VPK087	Bowl Kit (Short)							No OEM Part # Available
VPK088	Bowl Kit (Long)							No OEM Part # Available
VPN092	Clean Out-Nipple							No OEM Part # Available
VPS076	Screen (50 Mesh)	No OEM Part # Available	No OEM Part # Available	No OEM Part # Available	No OEM Part # Available	No OEM Part # Available	No OEM Part # Available	
VPS077	Screen (20 Mesh)	No OEM Part # Available	No OEM Part # Available	No OEM Part # Available	No OEM Part # Available	No OEM Part # Available	No OEM Part # Available	
VPS078	Screen (50 Mesh)	No OEM Part # Available	No OEM Part # Available	No OEM Part # Available	No OEM Part # Available	No OEM Part # Available		
VPS079	Screen (20 Mesh)	No OEM Part # Available	No OEM Part # Available	No OEM Part # Available	77005079	77005079		
VPS089	Screen (Short)							No OEM Part # Available
VPS090	Screen (Long)		DIALEOGEO					No OEM Part # Available
VPV049	Check Valve (1/2")	56105	PW50656	64569101	770000000	77000000		
VPV051 VPV053	Check Valve (3/4") Check Valve (1")	56105	PVV50657 PVV50658	64568191 64568193	770002003	770002003		
VPV053 VPV055	Check Valve (1) Check Valve (3/8")		1 1 1 1 1 1 1 1 1 1 1 1 1 1 1 1 1 1 1 1	64622311				
VPV055 VPV057	Check Valve (1/4" In-Line)	56056		07022311				
VPV057 VPV058	Check Valve (3/4")	00000					CVF-750	
VPV059	Vacuum Relief Valve	55865						
VPV060	Vacuum Relief Valve (1")		PW85980 & PW85981					
VPV063	Vacuum Relief Valve (3/4")		F V V O D D O I				VRV-100	
VPV063 VPV083	Cuspidor Drain Valve (3/4)						VIIV-100	4512266
	caopiaol Brain valvo (vala hap)	I		I	1		1	4012200