




































FIXED PATIENT CABLES

RPI PART #	DESCRIPTION	# OF LEADS	LENGTH	EQUIPMENT END	PATIENT END	NOTES
KCB002	PATIENT CABLE - 3 LEAD FIXED SNAP	3	~11 ft.	6-Pin; 1K Ω Resistors Jumper E to C 	Fixed Snap 	
KCB008	PATIENT CABLE - 3 LEAD FIXED SNAP	3	~11 ft.	6-Pin Jumper E to B 	Fixed Snap 	
KCB036	PATIENT CABLE - 3 LEAD FIXED SNAP	3	~10 ft.	12-Pin Right Angle 	Fixed Snap 	<ul style="list-style-type: none"> • OEM Part #11110-00029 & 3006218-02
KEB012	PATIENT CABLE - 5 LEAD FIXED SNAP	5	~13 ft.	6-Pin 	Fixed Snap 	
KEB014	PATIENT CABLE - 5 LEAD FIXED PIN	5	~8 ft.	5-Pin 	Fixed .125" Pin 	<ul style="list-style-type: none"> • This cable replaces RPI Part #BUC009 • Also available, Leadwire Adapter Kit (RPI Part #LAK001)
KEB030	PATIENT CABLE - 5 LEAD FIXED BANANA	5	~13 ft.	6-Pin 	Fixed Banana Plug 	<ul style="list-style-type: none"> • This cable replaces RPI Part #HPC005 • Also available, Leadwire Adapter Kit (RPI Part #LAK001)
KKB002	PATIENT CABLE - 10 LEAD FIXED BANANA	10	~10 ft.	15-Pin D-Sub 	Fixed Banana Plug 	<ul style="list-style-type: none"> • Also available, Leadwire Adapter Kit (RPI Part #LAK001)

FIXED PATIENT CABLES *(Continued)*

RPI PART #	DESCRIPTION	# OF LEADS	LENGTH	EQUIPMENT END	PATIENT END	NOTES
KKB003	PATIENT CABLE - 10 LEAD FIXED BANANA	10	~10 ft.	15-Pin D-Sub 	Fixed Banana Plug  GREEN RED 	<ul style="list-style-type: none"> • Also available, Leadwire Adapter Kit (RPI Part #LAK001)
KKB006	PATIENT CABLE - 10 LEAD FIXED BANANA	10	~10 ft.	15-Pin D-Sub; 20K Ω Resistors 	Fixed Banana Plug  GREEN RED 	<ul style="list-style-type: none"> • Also available, Leadwire Adapter Kit (RPI Part #LAK001)
KKB007	PATIENT CABLE - 10 LEAD FIXED BANANA	10	~10 ft.	15-Pin D-Sub; 10K Ω Resistors 	Fixed Banana Plug  GREEN RED 	<ul style="list-style-type: none"> • Also available, Leadwire Adapter Kit (RPI Part #LAK001)
KKB008	PATIENT CABLE - 10 LEAD FIXED BANANA	10	~10 ft.	15-Pin D-Sub; 10K Ω Resistors 	Fixed Banana Plug  GREEN RED 	<ul style="list-style-type: none"> • Also available, Leadwire Adapter Kit (RPI Part #LAK001)
KKB010	PATIENT CABLE - 10 LEAD FIXED PINCH	10	~18 ft.	12-Pin; 10K Ω Resistors 	Fixed Pinch  GREEN RED 	<ul style="list-style-type: none"> • OEM Part #60-00180-01 & 60-00184-01
KKB011	PATIENT CABLE - 10 LEAD FIXED SNAP	10	~18 ft.	12-Pin ; 10K Ω Resistors 	Fixed Snap  GREEN RED 	<ul style="list-style-type: none"> • OEM Part #60-00181-01 & 60-00185-01
KKB012	PATIENT CABLE - 10 LEAD FIXED BANANA	10	~9 ft.	15-Pin D-Sub; 10K Ω Resistors 	Fixed Banana Plug  GREEN RED 	<ul style="list-style-type: none"> • OEM Part #223 418 09 • Also available, Leadwire Adapter Kit (RPI Part #LAK001)