




## HOLTER CABLES

RPI PART #	DESCRIPTION	# OF LEADS	LENGTH	EQUIPMENT END	PATIENT END	SUGGESTED LEADWIRE KITS
<b>KEA004</b>	HOLTER CABLE - 5 LEAD 3 CHANNEL	5	~30"	7-Socket Jumper 1 to 5 and 2 to 6	Din/Unshielded .060" Pin (Safety) 	<b>Reusable Leadwire Sets</b> RPI Part #'s LKM026, LKM029 & LKM030 <b>Disposable Leadwire Sets</b> RPI Part #'s LDM303 & LDM310 <b>Veterinary Leadwire Set</b> RPI Part #LVM703
<b>KEB020</b>	HOLTER CABLE - 5 LEAD FIXED SNAP	5	~43"	7-Socket	Fixed Snap 	<b>NOTES:</b> • OEM Part #3-370-0001
<b>KGA001</b>	HOLTER CABLE - 7 LEAD 3 CHANNEL	7	~40"	7-Socket	Din/Unshielded .060" Pin (Safety) 	<b>Reusable Leadwire Sets</b> RPI Part #'s LKM038 & LKM041 <b>Disposable Leadwire Sets</b> RPI Part #'s LDM304 & LDM311 <b>Veterinary Leadwire Set</b> RPI Part #LVM704