




RPI is your Source for parts to fit Pelton & Crane Magnaclave Sterilizers!

PELTON & CRANE



Kit Includes all parts as shown. All parts also sold separately!

RPI Part #PCK240
OEM Part #004598, 004864, 004599, 005929, & 091007


DRAIN KIT
Fits: Drain Adapter and Tube Assembly

Parts also sold separately:
RPI Part #PCC241
OEM Part #004598 (Cap) & 004864 (O-Ring)
DRAIN CAP WITH O-RING
• Includes O-Ring (RPI Part #RPO386)

RPI Part #RPO386 (OEM Part #004864)
O-RING - 6/PKG

RPI Part #RPH285 (No OEM Part # Available)
DRAIN HOSE - Sold by the foot

RPI Part #RPC442 (No OEM Part # Available)
UNICLAMP - 6/PKG



Kit Includes all parts as shown. All parts also sold separately!


RPI Part #PCK238
OEM Part #004630 & A65-004630
CONDENSATION COIL KIT
Fits: Reservoir

Parts also sold separately:
RPI Part #RPF221 (No OEM Part # Available)
COMPRESSION NUT - 5/pkg

RPI Part #RPF217 (No OEM Part # Available)
COMPRESSION SLEEVE - 5/pkg

RPI Part #PCK226 (No OEM Part # Available)
FILTER WITH COTTER PIN


The Filter With Cotter Pin (RPI Part #PCK226) is also included in the Sterilizer PM Kit (RPI Part #PCK216) – see page 2.




RPI Part #PCK236
OEM Part #5473041
SOLENOID KIT (DOOR LATCH)
Fits: Door Locking Bracket

Part also sold separately:
RPI Part #PCS239
OEM Part #004525 & A63-004525
SPRING (DOOR INTERLOCK)
Fits: Door Lock Pin

RPI Part #PCK242
OEM Part # (No OEM Part # Available)
BUSHING KIT (DOOR HINGE)
Fits: Door Support



Includes all parts as shown!




RPI Part #PCB237
OEM Part #004800 (Bellows only)

BELLOWS (PRESSURE SWITCH)
Fits: Pressure Switch and Frame Assembly

RPI Part #PCS094
OEM Part #004478

PRESSURE SWITCH
• Mounts on frame of bellows assembly




RPI Part #PCS235
OEM Part #004837 & B02-004837

RELAY SOCKET
• 8 Pin
Fits: Reservoir Support

RPI Part #PCR196
OEM Part #004463

LATCHING RELAY
• 8 pin octal configuration
Fits: Relay Socket



RPI Part #PCR234
OEM Part #041002 & A62-011081

RELAY (ELECT. CONTACTOR)
Fits: Rear Chassis Support



replacement parts industries, inc.


CALL (800) 221-9723 • FAX (818) 882-7028 • EMAIL order@rpiparts.com • WEBSITE www.rpiparts.com

The above parts are manufactured by Replacement Parts Industries, Inc. to fit Pelton & Crane equipment. All product names used in this document are trademarks or registered trademarks of their respective holders. Web 03/11

Table Top Sterilizers



Parts to fit Pelton & Crane Magnaclave Sterilizers



Kit Includes all parts as shown.
All parts also sold separately!

RPI Part #PCK216
OEM Part #004497 (Door Gasket Only)
STERILIZER PM KIT

Parts also sold separately:
RPI Part #PCG019 (OEM Part #004497)
DOOR GASKET
• 16.500" OD
Fits: Door

RPI Part #PCB001 (OEM Part #004048)
BELLOWS
Fits: Air Valve Housing


RPI Part #PCK226 (No OEM Part # Available)
FILTER WITH COTTER PIN
Fits: Condensation Coil

RPI Part #PCF009 (OEM Part #1881031 & 004326)
FILTER (FILL LINE)
Fits: Fill Tube

RPI Part #RPO354 (OEM Part #004004)
AIR VALVE O-RING - 6/PKG
• .926" ID x 1.066" OD x .070" C/S; Silicone
Fits: Air Valve Housing Bellows Cap



Download your Free
Sterilizer PM Checklist
from the RPI Tech Help page on
our website – www.rpiparts.com

RPI Part #PCT013
OEM Part #004640
FILL VENT TUBE
Fits: Fill Filter and Chamber



RPI Part #RPT113
OEM Part # (No OEM Part # Available)

MAX REGISTER THERMOMETER
• Mercury; Length: 5-1/2"
• Use for temperature verification in sterilizer




Kit Includes all parts as shown.
All parts also sold separately!

RPI Part #PCK224
OEM Part #011047
BELLOWS KIT

Parts also sold separately:
RPI Part #PCB001 (OEM Part #004048)
BELLOWS
Fits: Air Valve Housing

RPI Part #RPO354 (OEM Part #004004)
AIR VALVE O-RING - 6/PKG
• .926" ID x 1.066" OD x .070" C/S; Silicone
Fits: Air Valve Housing Bellows Cap


RPI Part #PCC005 (OEM Part #004228)
BELLOWS CAP
• Requires Air Valve O-Ring (RPI Part #RPO354)
Fits: Air Valve Housing



Includes all parts as shown!

RPI Part #PCG040
OEM Part #3336356 & 015316


PRESSURE GAUGE
• Dual Scale: 0 to 60 psi / 0 to 400 kpa
Fits: Control Panel



Includes all parts as shown!

RPI Part #PCG079
OEM Part #4207796 & 015363

TEMPERATURE GAUGE
• Dual scale: 100°F to 330°F / 94°C to 165°C
• Capillary tubing length: 54"
Fits: Control Panel



Includes Terminal Hardware!

RPI Part #PCH091
OEM Part #004515

MAIN HEATER STRIP (208VAC)
• 1 KW, 208 VAC, ~43Ω
Fits: Chamber

RPI Part #PCH089
OEM Part #004516

WALL HEATER STRIP (208VAC)
• 500W, 208 VAC, ~87Ω
Fits: Chamber

RPI Part #PCH083
OEM Part #004502


MAIN HEATER STRIP (240VAC)
• 1 KW, 240 VAC, ~58Ω
Fits: Chamber

RPI Part #PCH092
OEM Part #004503

WALL HEATER STRIP (240VAC)
• 500 W, 240 VAC, ~115Ω
Fits: Chamber


RPI Part #PCT197
OEM Part #004505

OVERHEAT THERMOSTAT SWITCH
• Switch opens @ 470°F ±20°F
Fits: Mounting hole on heating element




RPI Part #PCT082
OEM Part #004527

WALL THERMOSTAT
• Opens @ 245° F,
Closes @ 215° F



RPI Part #PCV109
OEM Part #004529

PRESSURE RELIEF VALVE
Fits: Chamber




Parts to fit Pelton & Crane Magnaclave Sterilizers




RPI Part #PCT081
OEM Part #004531 & 041001

ELECTRIC TIMER (60 MINUTE)
Fits: Control Panel




RPI Part #PCB088
OEM Part #004491

BUZZER
Fits: Mounts on machine Base




RPI Part #PCS043
OEM Part #004079

FUNCTION SWITCH (3 LEAD)
• Mounts to bracket in main valve assembly




RPI Part #PCS044
OEM Part #004073

POWER SWITCH (2 LEAD)
• Mounts to bracket in main valve assembly




RPI Part #PCS090
OEM Part #004545

DOOR INTERLOCK SWITCH
• Mounts to bracket on door locking mechanism




RPI Part #PCS086
OEM Part #004468

DOOR OPEN SWITCH
• Mounts to bracket on partition in reservoir area



RPI Part #PCS087
OEM Part #004461

1 LB. PRESSURE SWITCH
• Mounts to bracket on one pound pressure switch assembly




RPI Part #PCN136
OEM Part #004076

SPEEDNUT - 10/PKG
• "U" type #10-24 thread



Valve supplied with "Packings" pressed in!



Includes all parts as shown!

RPI Part #PCV075
OEM Part #004310

MAIN VALVE BODY
• Valve supplied with "packings" pressed-in
Fits: Main Valve Assembly

RPI Part #PCC073
OEM Part #004030

MAIN VALVE CAP
Fits: Main Valve




RPI Part #PCN072
OEM Part #004027

MAIN VALVE NUT
Fits: Main Valve



RPI Part #PCS074
OEM Part #004028

MAIN VALVE STEM
Fits: Main Valve



TO REPAIR VALVE,
YOU WILL NEED
(2) KITS!



Kit Includes all parts as shown.
Parts also sold separately!

RPI Part #PCK026
OEM Part #017706 (OEM has only 1 set)


MAIN VALVE KIT
• To repair valve, need (2) Main Valve Kits
(RPI Parts #PCK026)
Fits: Main Valve

Parts also sold separately:
RPI Part #PCR037 (OEM Part #004039)
VALVE BUTTON RETAINER
Fits: Main Valve

RPI Part #RPO355 (OEM Part #004000)
MAIN VALVE O-RING - 6/PKG
• .801" ID x .941" OD x .070" C/S; Silicone
Fits: Main Valve


RPI Part #PCA118
OEM Part #004231

LOCATOR ARM - 3/pkg
• Mounts to Main Valve




RPI Part #PCL029
OEM Part #004501

LAMP (OPEN DOOR)
Fits: Control Panel




RPI Part #PCL030
OEM Part #004499

LAMP (HEAT ON)
Fits: Control Panel




RPI Part #PCL031
OEM Part #004500

LAMP (STERILE)
Fits: Control Panel



RPI Part #PCL032
OEM Part #004498

LAMP (POWER ON)
Fits: Control Panel





List of Parts to fit Pelton & Crane Magnaclave Sterilizers

Table with 3 columns: RPI Part #, OEM Part #, Description. Lists various parts like Locator Arm, Bellows, Buzzer, etc.

Table with 3 columns: RPI Part #, OEM Part #, Description. Lists various parts like Main Valve Nut, Speednut, Valve Button Retainer, etc.

Visit our website www.rpiparts.com for more information about these parts.

RPI is your Best Source for parts to fit Pelton & Crane Magnaclave sterilizers!

ADDRESS SERVICES REQUESTED

PRESORTED STANDARD U.S. POSTAGE PAID VAN NUYS, CA PERMIT NO. 2399

replacement parts industries, inc. 20338 Corisco Street • Chatsworth, California • 91311

