

RPI is Your Best Source for Table Top Sterilizer Parts!



- 28 different PM Kits for table top sterilizers serviced most
- Outstanding Technical Service – FREE via phone, fax and e-mail
- Troubleshooting Guides and PM Checklists – FREE downloads from www.rpiparts.com

PM Kit to fit AMSCO Eagle® 10 & 10 Plus




RPI Part #AMK234
OEM Part #
P426637-261
(Gasket Only)

STERILIZER PM KIT

- Includes Door Gasket (RPI Part #AMG161), Water Guard (RPI Part #AMG172), Filter (Fill Line) (RPI Part #PCF009), and a PM Check List

Models: Eagle 10 and Eagle 10 plus



PM Kits to fit Midmark® • Ritter® 7, M7, M9 & M11




RPI Part #MIK134
OEM Part #
002-0243-00
(Gasket Only)

STERILIZER PM KIT

- Includes Fill/Vent Mesh Chamber Filter (RPI Part #MIF062), Door Gasket (RPI Part #RCG019), and a PM Check List

Models: M7 and 7




PM Kits to fit SCICAN STATIM® 2000 & 5000

RPI Part #SCK035
OEM Part #
01-100028S
(Gasket Only)

STERILIZER PM KIT

- Includes Cassette Seal (RPI Part #SCS001), Liquid Soap (RPI Part #RPS287), Compressor Filter (RPI Part #SCF002), Biological Filter (RPI Part #SCF032 for S/N 2101DK0274 and above), Water Filter (In-Line) (RPI Part #SCF034), and a PM Check List

Model: 2000




RPI Part #MIK072
OEM Part #
002-0361-00
(Door & Dam Gaskets only)

STERILIZER PM KIT

- Includes Mesh Style Filters (RPI Part #MIF062), Sintered Bronze Filter (RPI Part #MIF038), Door Springs (RPI Part #MIS076), Door Gasket (RPI Part #MIG027), Dam Gasket (RPI Part #MIG035), and a PM Check List

Models: M9 & M9D




RPI Part #SCK036
OEM Part #
01-101649S
(Gasket Only)

STERILIZER PM KIT

- Includes Cassette Seal (RPI Part #SCS029), Liquid Soap (RPI Part #RPS287), Air Compressor Filter (RPI Part #SCF031 for use with MEDO Compressor only), Biological Filter (RPI Part #SCF032), Water Filter (In-Line) (RPI Part #SCF034), and a PM Check List

Model: 5000




RPI Part #MIK080
OEM Part #
002-0504-00
(Door & Dam Gaskets only)

STERILIZER PM KIT

- Includes Mesh Style Filters (RPI Part #MIF062), Sintered Bronze Filter (RPI Part #MIF038), Door Springs (RPI Part #MIS076), Door Gasket (RPI Part #MIG028), Dam Gasket (RPI Part #MIG036), and a PM Check List

Models: M11 & M11D



replacement parts industries, inc.

CALL (800) 221-9723 • FAX (818) 882-7028 • EMAIL order@rpiparts.com • WEBSITE www.rpiparts.com

The above parts are manufactured by Replacement Parts Industries, Inc. to fit the equipment noted. All product names used in this document are trademarks or registered trademarks of their respective holders. Web 01/11

STERILIZER PM KITS – AMERICAN STERILIZER, MIDMARK® • RITTER®, PELTON & CRANE, SCICAN, & TUTTNAUER®

SECTION

10

PAGE


116

UPDATED

01/11



PM Kits to fit Tuttnauer EZ9, EZ10, EZ10k, 1730, 2340, 2540 & 3870 series




RPI Part #TUK121
OEM Part #02610020
(Door Gasket Only)

STERILIZER PM KIT

- Includes Fill/Vent Mesh Chamber Filter (RPI Part #MIF062), Door Bellows (RPI Part #TUB019), Rubber Boot For Door Switch (RPI Part #TUB072), Door Gasket (RPI Part #TUG001), and Air Jet Valve (Black Top) (RPI Part #TUJ034), and a PM Check List

Models: 1730E/M




RPI Part #TUK126
OEM Part #02610118
(Door Gasket Only)

STERILIZER PM KIT

- Includes Fill/Vent Mesh Chamber Filter (RPI Part #MIF062), Door Bellows (RPI Part #TUB019), Rubber Boot For Door Switch (RPI Part #TUB072), Door Gasket (RPI Part #TUG021), Air Jet Valve (Red Top) (RPI Part #TUJ033), and a PM Check List

Models: 2340EK/MK




RPI Part #TUK122
OEM Part #02610020
(Door Gasket Only)

STERILIZER PM KIT

- Includes Fill/Vent Mesh Chamber Filter (RPI Part #MIF062), Door Bellows (RPI Part #TUB019), Rubber Boot For Door Switch (RPI Part #TUB072), Door Gasket (RPI Part #TUG001), and Air Jet Valve (Red Top) (RPI Part #TUJ033), and a PM Check List

Models: 1730EK/MK/MKV




RPI Part #TUK128
OEM Part #02610023
(Door Gasket Only)

STERILIZER PM KIT

- Includes Fill/Vent Mesh Chamber Filter (RPI Part #MIF062), Door Bellows (RPI Part #TUB019), Rubber Boot For Door Switch (RPI Part #TUB072), Door Gasket (RPI Part #TUG003), Air Jet Valve (Black Top) (RPI Part #TUJ034), and a PM Check List

Models: 2540E/M




RPI Part #TUK123
OEM Part #02610005
(Door Gasket Only)

STERILIZER PM KIT

- Includes Fill/Vent Mesh Chamber Filter (RPI Part #MIF062), Door Bellows (RPI Part #TUB019), Rubber Boot For Door Switch (RPI Part #TUB072), Door Gasket (RPI Part #TUG002), Air Jet Valve (Black Top) (RPI Part #TUJ034), and a PM Check List

Model: 2340M (Serial #8805 and below).




RPI Part #TUK130
OEM Part #02610023
(Door Gasket Only)

STERILIZER PM KIT

- Includes Fill/Vent Mesh Chamber Filter (RPI Part #MIF062), Door Bellows (RPI Part #TUB019), Rubber Boot For Door Switch (RPI Part #TUB072), Door Gasket (RPI Part #TUG003), Air Jet Valve (Red Top) (RPI Part #TUJ033), and a PM Check List

Models: 2540EK/MK




RPI Part #TUK124
OEM Part #02610118
(Door Gasket Only)

STERILIZER PM KIT

- Includes Fill/Vent Mesh Chamber Filter (RPI Part #MIF062), Door Bellows (RPI Part #TUB019), Rubber Boot For Door Switch (RPI Part #TUB072), Door Gasket (RPI Part #TUG021), Air Jet Valve (Black Top) (RPI Part #TUJ034), and a PM Check List

Models: 2340E & 2340M (Serial #8806 and above).




RPI Part #TUK131
OEM Part #02610019
(Door Gasket Only)

STERILIZER PM KIT

- Includes Fill/Vent Mesh Chamber Filter (RPI Part #MIF062), Door Bellows (RPI Part #TUB019), Rubber Boot For Door Switch (RPI Part #TUB072), Door Gasket (RPI Part #TUG074), Air Jet Valve (Black Top) (RPI Part #TUJ034), and a PM Check List

Models: 3870E/M

PM Kits to fit Tuttnauer EZ9, EZ10, EZ10k, 1730, 2340, 2540 & 3870 series - continued




RPI Part #TUK125
OEM Part #02610118
 (Door Gasket Only)

STERILIZER PM KIT

- Includes Fill/Vent Mesh Chamber Filter (RPI Part #MIF062), Door Bellows (RPI Part #TUB019), Rubber Boot For Door Switch (RPI Part #TUB072), Door Gasket (RPI Part #TUG021), and Air Jet Valve (Black Top) (RPI Part #TUJ034), HEPA Filter (RPI Part #RPF644), and a PM Check List

Models: EZ9 & 2340EA




RPI Part #TUK129
OEM Part #02610023
 (Door Gasket Only)

STERILIZER PM KIT

- Includes Fill/Vent Mesh Chamber Filter (RPI Part #MIF062), Door Bellows (RPI Part #TUB019), Rubber Boot For Door Switch (RPI Part #TUB072), Door Gasket (RPI Part #TUG003), Air Jet Valve (Red Top) (RPI Part #TUJ033), HEPA Filter (RPI Part #RPF644), and a PM Check List

Models: EZ10, EZ10k, 2540EA/EKA/EHS




RPI Part #TUK127
OEM Part #02610118
 (Door Gasket Only)

STERILIZER PM KIT

- Includes Fill/Vent Mesh Chamber Filter (RPI Part #MIF062), Door Bellows (RPI Part #TUB019), Rubber Boot For Door Switch (RPI Part #TUB072), Door Gasket (RPI Part #TUG021), Air Jet Valve (Red Top) (RPI Part #TUJ033), HEPA Filter (RPI Part #RPF644), and a PM Check List

Model: 2340EKA



RPI Part #TUK132
OEM Part #02610019
 (Door Gasket Only)


STERILIZER PM KIT

- Includes Fill/Vent Mesh Chamber Filter (RPI Part #MIF062), Door Bellows (RPI Part #TUB019), Rubber Boot For Door Switch (RPI Part #TUB072), Door Gasket (RPI Part #TUG074), Air Jet Valve (Black Top) (RPI Part #TUJ034), HEPA Filter (RPI Part #RPF644), and a PM Check List

Models: 3870EA/EHS

PM Kits to fit Pelton & Crane Delta AE8, XL8, Q8, AF10, XL10 & Q10

RPI Part #PCK231
OEM Part #1539241
 (Door Gasket Only)




STERILIZER PM KIT

- Includes Compressor Filter (RPI Part #SCF002), Door Gasket (RPI Part #PCG065), Filter (Fill Line) (RPI Part #PCF009), Filter In-Line Assembly (RPI Part #PCF233), HEPA Filter (RPI Part #RPF644), Reservoir Filter Kit (RPI Part #PCK167), Spacer (Door) (RPI Part #PCS067), and a PM Check List

Models: AE8, XL8 & Q8

RPI Part #PCK232
OEM Part #1539407
 (Door Gasket Only)



STERILIZER PM KIT


- Includes Compressor Filter (RPI Part #SCF002), Door Gasket (RPI Part #PCG066), Filter (Fill Line) (RPI Part #PCF009), Filter In-Line Assembly (RPI Part #PCF233), HEPA Filter (RPI Part #RPF644), Reservoir Filter Kit (RPI Part #PCK167), Spacer (Door) (RPI Part #PCS067), and a PM Check List

Models: AF10, XL10 & Q10

More PM Kits to fit Pelton & Crane on page 4!

Table Top Sterilizers


PM Kits to fit Pelton & Crane Magnaclave, OCM, OCR, OCR+, Sentry & Validator



RPI Part #PCK216
OEM Part #004497
 (Door Gasket Only)

STERILIZER PM KIT


- Includes Door Gasket (RPI Part #PCG019), Air Valve O-ring (RPI Part #RPO354), Bellows (RPI Part #PCB001), Filter with Cotter Pin (RPI Part #PCK226), Filter (Fill Line) (RPI Part #PCF009), and a PM Check List
- Model:** Magnaclave



RPI Part #PCK220
OEM Part #1539241
 (Door Gasket Only)

STERILIZER PM KIT


- Includes Door Gasket (RPI Part #PCG065), Filter (Fill Line) (RPI Part #PCF009), Air Valve O-ring (RPI Part #RPO354), Bellows (RPI Part #PCB001), Spacer (Door) (RPI Part #PCS067), and a PM Check List
- Model:** Validator 8



RPI Part #PCK217
OEM Part #001187
 (Door Gasket Only)

STERILIZER PM KIT


- Includes Door Gasket (RPI Part #PCG020), Filter (Fill Line) (RPI Part #PCF009), Air Valve O-ring (RPI Part #RPO354), Bellows (RPI Part #PCB001), and a PM Check List
- Model:** Sentry



RPI Part #PCK221
OEM Part #1539407
 (Door Gasket Only)

STERILIZER PM KIT


- Includes Door Gasket (RPI Part #PCG066), Filter (Fill Line) (RPI Part #PCF009), Air Valve O-ring (RPI Part #RPO354), Bellows (RPI Part #PCB001), Spacer (Door) (RPI Part #PCS067), and a PM Check List
- Model:** Validator 10



RPI Part #PCK218
OEM Part #004014
 (Door Gasket Only)

STERILIZER PM KIT


- Includes Door Gasket (RPI Part #PCG014), Filter (Fill Line) (RPI Part #PCF009), Air Valve O-ring (RPI Part #RPO354), Bellows (RPI Part #PCB001), and a PM Check List
- Model:** OCM



RPI Part #PCK222
OEM Part #1539241
 (Door Gasket Only)

STERILIZER PM KIT


- Includes Door Gasket (RPI Part #PCG065), Filter (Fill Line) (RPI Part #PCF009), Spacer (Door) (RPI Part #PCS067), Reservoir Filter Kit (RPI Part #PCK167), and a PM Check List
- Models:** Validator Plus 8



RPI Part #PCK219
OEM Part #004341
 (Door Gasket Only)

STERILIZER PM KIT

- Includes Door Gasket (RPI Part #PCG015), Filter (Fill Line) (RPI Part #PCF009), Air Valve O-ring (RPI Part #RPO354), Bellows (RPI Part #PCB001), and a PM Check List
- Model:** OCR & OCR+



RPI Part #PCK223
OEM Part #1539407
 (Door Gasket Only)

STERILIZER PM KIT

- Includes Door Gasket (RPI Part #PCG066), Filter (Fill Line) (RPI Part #PCF009), Spacer (Door) (RPI Part #PCS067), Reservoir Filter Kit (RPI Part #PCK167), and a PM Check List
- Models:** Validator Plus 10