




New Parts to fit


MIDMARK® • RITTER®
UltraClave® M9 & M11 Sterilizers!




RPI Part #MIH238
OEM Part #015-1554-00 (Main Harness)/015-1640-00 (Hi Limit Jumper)

MAIN HARNESS

- Includes: Both the main harness and the high limit jumper
- Fits:** Main Control Board
- Models:** M9-020 thru -022 (Serial # prefixes RN, RP, RR & V), M9D-020 thru -022 (Serial # prefixes RW, RX & V), M11-020 thru -022 (Serial # prefixes RS, RT, RV & V) & M11D-020 thru -022 (Serial # prefixes RY, RZ & V)




RPI Part #MIA180
OEM Part #002-0790-00
DOOR MOTOR ASSEMBLY
 • 120 VAC
Fits: Chamber Frame
Models: M9-020 thru 022 & M11-020 thru 022





RPI Part #MIC226
OEM Part #053-0505-00
PRINTER CARTRIDGE
 • Epson® brand ribbon cartridge; Color: Black
Fits: Printers (Old Style and New Style)
Models: M9-001 thru 019, M9-020 thru 022, M11-001 thru 019 & M11-020 thru 022



RPI Part #MIK236
OEM Part #002-1244-01
STEAM BLOCK KIT
 • Material: Beige silicon rubber
 • Size: 13.75" lg.
 • Includes all parts as shown
Fits: Door Panel
Models: M9-020 thru -022 & M9D-020 thru -022



MIK236
 Cross Section
 Shown Actual Size





RPC957
 (4 pcs)






RPI Part #RPC957
OEM Part #042-0182-00
SNAP CLIP
 • 10 per package
 • Fits Hole Dia.: .156", Panel Thk.: .093"
 • Material: Nylon; Color: Natural
 • Head Dia.: .312"
 • Included in: Steam Block Kit (RPI Part #MIK236 & MIK237)
Fits: Steam Block
Models: M9-020 thru -022, M9D-020 thru -022, M11-020 thru -022 & M11D-020 thru -022

RPI Part #MIK237
OEM Part #002-1244-02
STEAM BLOCK KIT
 • Material: Beige silicon rubber
 • Size: 16.25" lg.
 • Includes all parts as shown
Fits: Door Panel
Models: M11-020 thru -022 & M11D-020 thru -022



MIK237
 Cross Section
 Shown Actual Size

RPC957
 (4 pcs)




replacement parts industries, inc.

CALL (800) 221-9723 • FAX (818) 882-7028 • EMAIL order@rpiparts.com • WEBSITE www.rpiparts.com

The above parts are manufactured by Replacement Parts Industries, Inc. to fit Midmark®/Ritter® equipment. All product names used in this document are trademarks or registered trademarks of their respective holders. Flyer Code 1513MI Show 08/16


Don't Forget ... RPI Offers the Top Cover Kits to fit Midmark® UltraClave® M9 & M11 Sterilizers!



RPI Part #MIK196
OEM Part #002-0356-00

TOP COVER KIT
Models: M9-001 thru -019

Includes all parts as shown!



RPI Part #MIK197
OEM Part #002-0503-00

TOP COVER KIT
Models: M11-001 thru -019

Includes all parts as shown!

New Parts to fit Midmark® UltraClave® M9 & M11 Sterilizers!

ADDRESS SERVICES REQUESTED

PERMIT NO. 2399
VAN NUYS, CA
U.S. POSTAGE PAID
PRESORTED STANDARD

replacement parts industries, inc.
20338 Corisco Street • Chatsworth, California • 91311



Make Your Midmark® UltraClave® Look As Good As New with RPI's Top Cover Warning Labels!




RPI Part #MIL260
OEM Part #061-0309-00, 061-0309-01

LABEL (PR VALVE WARNING)

- Size: 1.50" x 4.21" with adhesive backing
- Included in: Top Cover Kits
(RPI Part #'s MIK196 & MIK197)

Fits: Top Cover Near the Inspection Cover
Models: M9-001 thru 019 & M11-001 thru 019



RPI Part #MIL261
OEM Part #061-0381-00, 061-0532-00

LABEL (WARNING)

- Size: 1.50" x 4.78" with adhesive backing
- Included in: Top Cover Kits
(RPI Part #'s MIK196 & MIK197)

Fits: Top Cover Near the Display Overlay
Models: M9-001 thru 019 & M11-001 thru 019

