



Wells Dental

Models: Engine Units: U900, U901, U903 & U904

Parts to fit Wells Dental!

Belts, Motor Brushes, Rubber Foot, Fuse and Power Cord


DENTAL LABORATORY EQUIPMENT

RPI Part #BDB001
OEM Part #E-107

BRUSH

- 2 per package
- Size: .250" x .187" x .660" lg.
- Spring Length: 1.10"
- Copper shunt wire

Models: Hood & Motor: E026, E027, E028 & E029 and Motor: Bodine (Small brush) - E-129 & E-130




RPI Part #WLB001
OEM Part #E-101

BRUSH

- 2 per package
- Size: .125" x .187" x .596" lg.
- Radius: .470"
- Spring Length: .82"
- Copper shunt wire

Models: Hood & Motor: E026, E027, E028 & E029 and Motor: Baldor - E132 & E133




RPI Part #WLB002
OEM Part #AS404

BELT (6' 6")

- Material: Pyrrathane
- Size: .093" dia. C/S x 6 feet & 6 inches circumference

Models: Arms: AS051 (Whisper short arm) and Arms: AS011 (Piccolo engine arm)




RPI Part #DRB001
OEM Part #E-105

BRUSH

- 2 per package
- Size: .250" x .375" x .817" lg.
- Radius: .625"
- Spring Length: 1.38"
- Copper shunt wire

Models: Hood & Motor: E026, E027, E028 & E029 and Motor: Bodine (Large brush) - E-131




RPI Part #EMB001
OEM Part #E-106

BRUSH

- 2 per package
- Size: .250" x .250" x .657" lg.
- Radius: .500"
- Spring Length: 1.34"
- Copper shunt wire

Models: Hood & Motor: E026, E027, E028 & E029 and Motor: WSM





RPI Part #WLB003
OEM Part #AS406

BELT (9' 8")

- Material: Pyrrathane
- Size: .093" dia. C/S x 9 feet & 8 inches circumference


Models: Arms: AS052 (Whisper long arm) and Arms: AS012 (Standard engine arm)



RPI Part #CAF011
OEM Part # (No OEM Part # Available)

RUBBER FOOT

- 4 per package
 - Material: Black rubber with molded in washer
 - Size: .750" dia. x .500" lg.
 - Mounts to base plate using #10 screw
- Models:** Hood & Motor: E026, E027, E028 & E029 and Foot Controls: F011 & F013




RPI Part #RPF051
OEM Part # (No OEM Part # Available)

FUSE (5A, 250V)

- 5 per package
- Type: Glass; AGC Fast Acting
- Size: 1/4" dia. x 1-1/4" lg.
- Also included in: Fuse Kit (RPI Part #RPK410)

Fits: Fuse Holder
Models: Foot Controls: F011 & F013



RPI Part #PCC004
OEM Part # (No OEM Part # Available)

POWER CORD WITHOUT CONNECTORS

- Length: 8 ft.; Color: Black
- 16/3 SJTW, 105°C; 13A@125VAC
- Plug Type: NEMA 5-15P
- Connector Type: None; .33" OD cord
- Use Strain Relief Bushing (RPI Part #RPB164) when replacing power cord
- Also available: Wire Nut (RPI Part #RPT017)

Models: Hood & Motor: E026, E027, E028 & E029



replacement parts industries, inc.

CALL (800) 221-9723 • FAX (818) 882-7028 • EMAIL order@rpiparts.com • WEBSITE www.rpiparts.com

The above parts are manufactured by Replacement Parts Industries, Inc. to fit Pelton & Crane equipment. All product names used in this document are trademarks or registered trademarks of their respective holders. Web 05/12