


O-RING KIT FOR SERVICING STERIS SYSTEM 1 STERILE PROCESSING SYSTEM

RPI Part #SSK050

- Kit includes 15-compartment storage case (RPI Part #RPC616), O-ring Lubricant (RPI Part #RPL615) and O-ring Removal Tools (RPI Part #RPK581)
- All parts are also sold separately



RPI PART #	STERIS PART #	X-SEC	ID	OD	INDUSTRY PART #	MATERIAL	DURO	KIT QTY	FITS
RPO307	450617	.070	.114	.254	AS568B-006	EPDM	70	12	Small In-Air Manifold
RPO354	450609	.070	.926	1.066	AS568B-021	Silicone	70	12	Small Q.D. Post in Drip Pan
RPO388	450610	.070	.426	.566	AS568B-013	Silicone	70	12	HP Small Droque Inside
RPO389	450608	.103	1.049	1.255	AS568B-121	Silicone	70	12	Kepner CK-4 Outer
RPO390	450614	.070	1.114	1.254	AS568B-024	Silicone	70	12	Kepner CK-4 Inner
RPO450	450621	.051	.547	.649	5-001	Silicone	70	12	CK-1, CK-2 & CK-3
RPO605	450600	.103	1.174	1.380	AS568B-123	Silicone	70	12	Large Q.D. Post under Drip Pan
RPO606	450601	.139	.859	1.137	AS568B-212	Silicone	70	12	Large - In Droque
RPO607	450605	.070	1.239	1.379	AS568B-026	Silicone	70	12	Large Q.D. Post in Drip Pan
RPO608	450606	.070	.489	.629	AS568B-014	Silicone	70	12	Air Manifold Check Valve
RPO609	450611	.070	.114	.254	AS568B-006	Silicone	70	12	Air-Water Valve
RPO610	450612	.103	.924	1.130	AS568B-119	Silicone	70	12	Small Q.D. Post under Drip Pan
RPO612	450620	3mm	80mm	86mm	M3X80	Silicone	60	3	Sterile Filter Housing Cap
RPO613	450622	.139	3.484	3.762	AS568B-238	Viton	70	3	Prefilter Housing